

INTRODUCING...

AEM® AT YOUR FINGERTIPS

AEM® ELECTROSURGICAL
INSTRUMENTS

**THE NEW DISPOSABLE
FIXED-TIP ELECTRODE...
now available
in hand switching**

ERGONOMIC HANDLE
WITH NO SLIP GRIP

EASY TWO-PIECE
ASSEMBLY

ELECTRODES
INTERCHANGEABLE IN
THE STERILE FIELD

ENCISION®

**THE NEW DISPOSABLE FIXED-TIP ELECTRODE...
now available in hand switching**



ELECTRODES INTERCHANGEABLE IN STERILE FIELD,
AVAILABLE IN SIX POPULAR STYLES

New from Encision®, the disposable Fixed-Tip Electrode, now available in hand switching. The dFTE offers an ergonomically designed handle with an easy reach rotation knob and locking feature for optimal electrode control.

- Sleek, ergonomic design for physician comfort
- Electrodes interchangeable in the sterile field
- No slip grip for more control
- Compatible with Encision's disposable sheath



NO SLIP COMFORT GRIP HANDLE

EASY TWO-PIECE ASSEMBLY



EASY REACH ROTATION KNOB WITH LOCKING FEATURE



UNIVERSAL ADAPTER FOR HAND SWITCHING



EASY ONE TIME ASSEMBLY OF THE UNIVERSAL ADAPTER ALLOWS OPERATION FOR HAND AND FOOT SWITCHING

COMPATIBLE WITH AEM® MONITORS ONLY

PRODUCT NO.	DESCRIPTION		
ES 1300	Hand Switching dFTE		
ES 9015	Universal Adapter		
dFTE ELECTRODES	DESCRIPTION	35 CM	45 CM
ES 0301	Spatula	ES 0301	ES 0301-45
ES 0309	Flat J-Hook	ES 0309	ES 0309-45
ES 0310	J-Hook	ES 0310	ES 0310-45
ES 0311	L-Diamond	ES 0311	ES 0311-45
ES 0312	L-Wedge	ES 0312	ES 0312-45
ES 0314	Square	ES 0314	ES 0314-45

Only Encision's instruments have AEM® TECHNOLOGY... designed to prevent stray energy burns caused by insulation failure or capacitive coupling in electrosurgery.

SAFE AND EFFECTIVE DELIVERY OF ELECTROSURGICAL ENERGY