



Reusable Modular Instruments

HOW TO ORDER



1 Ensure AEM® Monitor and ESU Compatibility.

- An AEM Monitor must be used in order to operate AEM monopolar laparoscopic instruments.

2 Select a Reusable Handle Style and a Shaft Length

- Select a handle style, either the ES8000 or the ES8200 Series enTouch® Handles.
- For a standard shaft length of 35cm, no suffix is added to the catalog number. For a shaft length of 45cm, the suffix "L" is added to the catalog number. For example, ES8000L.
- For the hooded handle with a standard shaft length of 35cm, the suffix 'H' is added to the catalog number. For example, ES8000H. For the hooded handle style with a shaft length of 45cm, the suffix 'HL' is added to the catalog number. For example, ES8000HL. The hooded handles are CE Marked for distribution in the EU and Australia only.

3 Select a Reusable Insert Style and a Shaft Length for the Reusable Insert

- Select a reusable insert style from the list of reusable graspers, dissectors, scissors, biopsy punches or spoons. For example, ES0013 Bullet Nose Grasper.
- For a standard shaft length of 35cm, no suffix is added to the catalog number. The catalog number would remain ES0013.
- For a shaft length of 45cm, the suffix -45 is added to the catalog number so the catalog number would become ES0013-45.

For Example:

Monitor		
Catalog #	EM2+	
Handle w/shaft	Length 35cm	Length 45cm
	Catalog #	ES8000 ES8000L
Insert	Length 35cm	Length 45cm (add '-45')
	Catalog #	ES0013 ES0013-45

Contact a sales representative for AEM Monitor and ESU compatibility, or to find an optimal insert pattern.

Precision Bipolar AEM® Monopolar

ENCISION




Reusable Articulating Inserts



PRODUCT HIGHLIGHTS

- Inserts must be used in conjunction with an ES8000 or an ES8200 enTouch® Handle and an AEM® Monitor.
- Inserts can be interchanged in the sterile field.
- Scissors Inserts can be used with the Disposable Sheath (ES0150A and ES0150-45).
- Distal tip configurations are available in two lengths (35cm or 45cm).



Reusable Scissors Inserts

Tip Styles	Jaw Length	Catalog #	
		35cm	45cm
  Curved Scissors, 1/2"	10mm	ES0001	ES0001-45
  Curved Scissors, 3/4"	17mm	ES0002	ES0002-45
  Hook Scissors	9mm Single Action	ES0010	ES0010-45

AEM® Monopolar

Bipolar

Precision

ENCISION

