


e•Access™ Slide Lock Handle

Accepts any Encision
Dissector or Grasper *


- Articulating rotation knob
- Locks when force is applied to the handle



Removable insert for ease
of cleaning

Integrated slide lock

1/4" thick hardened stainless
steel handles provide power



-
- Clean and sterilize with the insert installed or disassembled.
 - Entire lumen can be inspected for cleanliness.
 - Validated for manual and automated cleaning processes.
 - Validated for Gravity and Pre-Vac steam sterilization.
-

Cleaning brush can be inserted through entire length of handle



e•Access Slide Lock Handle



+

Reusable Insert



Order Information

		Length	Catalog #
e•Access™ Slide Lock Handle (rotatable)		35 cm	ES8500SLR
		45 cm	ES8500-45SLR
Insert, Atraumatic Grasper (Wave)		35 cm	ES0537
		45 cm	ES0537-45
Insert, Atraumatic Grasper (Single Action)		35 cm	ES0570
		45 cm	ES0570-45
Insert, Traumatic Grasper		35 cm	ES0571
		45 cm	ES0571-45
Insert, Traumatic Grasper (Single Action)		35 cm	ES0573
		45 cm	ES0573-45
Insert, Fenestrated Bowel Grasper		35 cm	ES0574
		45 cm	ES0574-45



Using a smart phone or tablet, the QR code provides quick and easy access to a setup video for the e•Access Handle.