

# **AEM® Disposable Electrodes Foot Control**

*The FC03XX Series Foot Control AEM Disposable Electrodes are single-piece, single-use, 5mm active instruments designed to fit commonly used 5mm trocar cannulae. AEM Disposable Electrodes are intended for ablation, removal, resection, and coagulation of soft tissue where associated hemostasis is required in open, endoscopic, and laparoscopic surgical procedures.*

## **Product Highlights**

- *AEM Disposable Electrodes are modular, single-use devices.*
- *Electrodes must be used in conjunction with an AEM Monitor.*
- *Electrodes can be used with a reusable or a disposable AEM Cord.*
- *Electrodes are shipped sterile and can be interchanged in the sterile field.*
- *Distal tip configurations are available in 35cm or 45cm lengths.*
- *Proximal cord connection rotates to accommodate easier target application.*
- *No handle assembly required; electrode connects directly to AEM Cord.*
- *Electrodes are available 10 per box.*

SAFE AND EFFECTIVE DELIVERY OF ELECTROSURGICAL ENERGY

# AEM® Disposable Electrodes, Foot Control

## How To Order



1


+

2



1. Ensure AEM Monitor and ESU compatibility.
  - An AEM Monitor must be used in order to operate AEM monopolar laparoscopic instruments.
2. Select an electrode tip style, for example, a Spatula Electrode is FC0301.
  - For a standard shaft length of 35cm, no suffix is added to the catalog number. The catalog number would remain FC0301.
  - For a shaft length of 45cm, the suffix -45 is added to the catalog number so the catalog number would become FC0301-45.

Note: A handle is not required for Foot Control AEM Disposable Electrodes. The AEM Cord connects directly to the electrode.

Tip Styles		Catalog #	
5mm (Box of 10)	Description	35cm	45cm
	AEM Disposable Electrode, Spatula, Foot Control	FC0301	FC0301-45
	AEM Disposable Electrode, Flat J-Hook, Foot Control	FC0309	FC0309-45
	AEM Disposable Electrode, J-Hook, Foot Control	FC0310	FC0310-45
	AEM Disposable Electrode, L-Diamond, Foot Control	FC0311	FC0311-45
	AEM Disposable Electrode, L-Wedge, Foot Control	FC0312	FC0312-45
	AEM Disposable Electrode, Square-L, Foot Control	FC0314	FC0314-45

SAFE AND EFFECTIVE DELIVERY OF ELECTROSURGICAL ENERGY