

AEM® Disposable Electrodes Foot Control

The FC03XX Series Foot Control AEM Disposable Electrodes are single-piece, single-use, 5mm active instruments designed to fit commonly used 5mm trocar cannulae. AEM Disposable Electrodes are intended for ablation, removal, resection, and coagulation of soft tissue where associated hemostasis is required in endoscopic and laparoscopic surgical procedures.



PRODUCT HIGHLIGHTS

- AEM Disposable Electrodes are modular, single-use devices.
- Electrodes must be used in conjunction with an AEM Burn Protection System.
- Electrodes are shipped sterile and can be interchanged in the sterile field.
- Distal tip configurations are available in 35cm or 45cm lengths.
- Proximal cord connection rotates to accommodate easier target application.
- No handle assembly required; electrode connects directly to AEM Cord.
- Electrodes are available 10 per box.

AEM® Disposable Electrodes Foot Control

HOW TO ORDER



5mm (Box of 10)	Description	35cm	45cm
	AEM Disposable Electrode, Spatula, Foot Control	FC0301	FC0301-45
	AEM Disposable Electrode, Flat J-Hook, Foot Control	FC0309	FC0309-45
	AEM Disposable Electrode, J-Hook, Foot Control	FC0310	FC0310-45
	AEM Disposable Electrode, L-Diamond, Foot Control	FC0311	FC0311-45
	AEM Disposable Electrode, L-Wedge, Foot Control	FC0312	FC0312-45
	AEM Disposable Electrode, Square-L, Foot Control	FC0314	FC0314-45

*The AEM EndoShield is currently undergoing premarket review by the FDA.