

## AEM® Disposable Electrodes Hand Control

The AEM Disposable Electrodes consist of the AEM Disposable Handswitching Handpiece with Rocker Switch and an integral cord and separate disposable electrodes that can be assembled in the sterile field. AEM Disposable Electrodes are intended for ablation, removal, resection, and coagulation of soft tissue where associated hemostasis is required in open, endoscopic, and laparoscopic surgical procedures.



### PRODUCT HIGHLIGHTS

- AEM Disposable Electrodes are modular, single use devices.
- The instrument cord is integral to the handpiece and does not require a separate AEM Cord.
- Electrodes can be used with Disposable Sheath (ES0155 and ES0155-45).
- Distal tip configurations are available in 35cm or 45cm lengths.
- Handpieces are available 25 per box.
- Electrodes are available 10 per box.

# AEM® Disposable Electrodes Hand Control

## AEM Disposable Handpiece

Tip Styles (Box of 25)	Catalog #
---------------------------	-----------



AEM Disposable  
Handswitching  
Handpiece with Rocker  
Switch

ES1300



## AEM Disposable Electrodes

Tip Styles 5mm (Box of 10)	Catalog #	
	35cm	45cm
 Spatula	ES0301	ES0301-45
 Flat J-Hook	ES0309	ES0309-45
 J-Hook	ES0310	ES0310-45
 L-Diamond	ES0311	ES0311-45
 L-Wedge	ES0312	ES0312-45
 Square-L	ES0314	ES0314-45