

## Instrument Accessories

### AEM® Cords

The reusable and disposable AEM® Cords are for use with AEM foot controlled instruments.


Cords		Catalog #
	AEM® Cord, Disposable	ES4107 (Box of 25)
	AEM® Cord, Reusable	ES4007

### Disposable Sheath

Disposable sheath is for use with reusable Scissors Inserts and reusable Fixed-Tip Electrodes.


- Fits over the Scissors Inserts and over the Fixed-Tip Electrode to control the area of exposed metal at distal end of instrument.
- Disposable sheaths available 10 per box.
- Disposable sheaths available for 35cm and 45cm length scissors inserts and electrodes.



Sheath		Outside Diameter	Catalog #	
			35cm	45cm
	Disposable Sheath (Box of 10)	5.5mm	ES0150A	ES0150-45

### Contact Cleaning Brush

- Cleaning brush is for use with is for use with AEM foot controlled instruments.
- It is designed to fit over the gold electrical contacts of an AEM instrument to help remove residue left on the contact during the normal cleaning and decontamination of the instrument.

Cleaning Brush		Catalog #
	Contact Cleaning Brush	BR1000

