

Reusable Fixed-Tip Electrodes

The ES3500 Series Fixed-Tip Electrodes are single-piece, 5mm active instruments designed to fit commonly-used 5mm trocar cannulae.



PRODUCT HIGHLIGHTS

- Electrodes must be used in conjunction with an AEM® Burn Protection System.
- Electrodes can be interchanged in the sterile field.
- Electrodes can be used with Disposable Sheath (ES0150A).
- Proximal cord connection rotates to accommodate easier target application.
- No handle assembly required; electrode connects directly to AEM® Cord.

Reusable Fixed-Tip Electrodes

Reusable Fixed-Tip Electrodes

Tip Styles		Catalog #
5mm		35cm
	Short Spatula	ES3501B
	Ball Tip	ES3504B
	Flat J-Hook	ES3509B
	J-Hook	ES3510B
	L-Diamond	ES3511B
	L-Wedge	ES3512B
	L-Hook	ES3513B
	Square-L	ES3514B
	Button Tip	ES3520B
	Needle Tip	ES3521B
	Corson Needle	ES3523B