


Reusable Modular Instruments

HOW TO ORDER



1 Ensure AEM® System and ESU Compatibility.

- An AEM Burn Protection System must be used in order to operate AEM monopolar laparoscopic instruments.

2 Select a Reusable Handle Style and a Shaft Length

- Select a handle style, either the ES8000 or the ES8200 Series enTouch® Handles.
- For a standard shaft length of 35cm, no suffix is added to the catalog number. For a shaft length of 45cm, the suffix "L" is added to the catalog number. For example, ES8000L.
- For the hooded handle with a standard shaft length of 35cm, the suffix 'H' is added to the catalog number. For example, ES8000H. For the hooded handle style with a shaft length of 45cm, the suffix 'HL' is added to the catalog number. For example, ES8000HL. The hooded handles are CE Marked for distribution in the EU and Australia only.

3 Select a Reusable Insert Style and a Shaft Length for the Reusable Insert

- Select a reusable insert style from the list of reusable graspers, dissectors, scissors, biopsy punches or spoons. For example, ES0013 Bullet Nose Grasper.
- For a standard shaft length of 35cm, no suffix is added to the catalog number. The catalog number would remain ES0013.
- For a shaft length of 45cm, the suffix -45 is added to the catalog number so the catalog number would become ES0013-45.

For Example:















Monitor		
Catalog #	EndoShield	
Handle w/shaft	Length 35cm	Length 45cm
Catalog #	ES8000	ES8000L
Insert	Length 35cm	Length 45cm (add '-45')
Catalog #	ES0013	ES0013-45

*The AEM EndoShield is currently undergoing premarket review by the FDA.


Reusable Dissector and Grasper Inserts

Tip Styles			Jaw Length	Catalog #	
				35cm	45cm
5mm (unless noted otherwise)					
		Curved Maryland Dissector	17mm	ES0501	ES0501-45
		Tapered Maryland Dissector	17mm	ES0526	ES0526-45
		Maryland Dissector, Aggressive	17mm	ES0549	ES0549-45
		Maryland Dissector, Diamond Serrations	17mm	ES0547	
		Maryland Dissector with 90° Tooth	19mm	ES0565	
		Beveled Maryland Dissector, Tapered	17mm	ES0566	ES0566-45
		Strong Curved Maryland, 7-8mm	19mm	ES0533	
		Strong Curved Maryland, 9-10mm	19mm	ES0531	
		Kelly Forceps	16mm	ES0541	ES0541-45











Reusable Dissector and Grasper Inserts

Tip Styles			Jaw Length	Catalog #	
				35cm	45cm
<i>5mm (unless noted otherwise)</i>					
		Right Angle Dissector	18mm	ES0004	
		Tapered Right Angle Dissector	19mm	ES0008	
		Dissecting Forceps, Right Angled	18mm	ES0548	
		Dissecting Forceps, Right Angled, Cross Serrated	18mm	ES0552	
		Short Right Angle Dissector	11mm	ES0011	
		Petelin Dissector	9mm Single Action	ES0510	
		Right Angle Dissector, 7mm, Diagonal Serrations	18mm	ES0559	



Reusable Dissector and Grasper Inserts

Tip Styles		Jaw Length	Catalog #	
5mm (unless noted otherwise)			35cm	45cm
		90° Grasper	20mm	ES0506
		Mixer Clamp, 90° Long	26mm	ES0538
		Straight Dissector	20mm	ES0513
		Pointed Nose Grasper	14mm	ES0509
		Bullet Nose Grasper	13mm	ES0013
		Bottle Nose Grasper	13mm	ES0012
		Dolphin Nose Grasper	17mm	ES0511
		Long Dolphin Nose Grasper	23mm	ES0512

Reusable Dissector and Grasper Inserts

Tip Styles			Jaw Length	Catalog #	
5mm (unless noted otherwise)				35cm	45cm
		Straight Micro-Grasper	11mm	ES0553	
		Duck Nose Micro-Grasper	9mm	ES0556	
		Curved Micro-Grasper	9mm	ES0554	
		Standard Grasper	13mm	ES0514	
		Blunt Nose Grasper	13mm	ES0009	ES0009-45
		Round Nose Grasper	12mm	ES0508	
		Fenestrated Grasping Forceps	19mm	ES0507	
		Fine Tooth Fenestrated Grasper	11mm	ES0014	

Reusable Dissector and Grasper Inserts

Tip Styles			Jaw Length	Catalog #	
				35cm	45cm
<i>5mm (unless noted otherwise)</i>					
		Micro-Fenestrated Grasper	10mm	ES0557	
		Tapered Micro-Fenestrated Grasper	10mm	ES0558	
		Endo Cinch Extreme Atraumatic Serrated	20mm	ES0535	
		Maxi Grasper	20mm	ES0543	ES0543-45
		Bowel Grasper	41mm	ES0521	ES0521-45
		Babcock-DeBaakey Forceps	31mm	ES0534	
		Ponce Forceps with Diamond Serrations	14mm	ES0560	
		Rodriguez Grasper	12mm Single Action	ES0561	

Reusable Dissector and Grasper Inserts

Tip Styles	Jaw Length	Catalog #	
5mm <i>(unless noted otherwise)</i>		35cm	45cm
	Uddo Grasper	17mm	ES0562

Reusable Biopsy Inserts

Tip Styles	Jaw Length	Catalog #	
5mm		35cm	45cm
	Biopsy Spoon	10mm	ES0551 ES0551-45