


## Bipolar Accessories



### End Point Monitor Remote Display

- For use with bipolar laparoscopic instruments only.
- Intended to aid the surgeon in determining the endpoint of bipolar electro-surgical desiccation through audible and visual communication.

Remote Display	Catalog #
 <p data-bbox="467 625 829 653">End Point Monitor Remote Display</p> <p data-bbox="496 667 764 695"><i>(for use with all AEM® Monitors)</i></p>	<p data-bbox="1040 646 1089 674">EMR</p>

### Reusable Bipolar Jumper Cords

- Required for use with End Point Monitoring.
- Used to connect the AEM® Monitor to the ESU.

Cords	Catalog #
 <p data-bbox="467 980 834 1008">Single Aspen Bipolar Jumper Cord</p> <p data-bbox="496 1022 976 1071"><i>(for use with EM2+ and 230V Covidien/Valleylab ESUs and selected 115V ConMed ESUs)</i></p>	<p data-bbox="1024 1014 1105 1041">ES9003</p>
 <p data-bbox="467 1121 792 1148">Standard Bipolar Jumper Cord</p> <p data-bbox="496 1163 964 1211"><i>(for use with EM2+ and selected 115V Covidien/Valleylab and ConMed ESUs)</i></p>	<p data-bbox="1024 1157 1105 1184">ES9004</p>