


e•Access™ Slide Lock Handle


PRODUCT HIGHLIGHTS

- Accepts any Encision Dissector or Grasper*
- Articulating rotation knob
- Locks when force is applied to the handle
- Removable insert for ease of cleaning
- 1/4" thick hardened stainless steel handles provide power

*Not compatible with other brands of inserts.



- Integrated slide lock



- Large flush port

- Large lumen for cleaning

- Clean and sterilize with the insert installed or disassembled.
- Entire lumen can be inspected for cleanliness.
- Validated for manual and automated cleaning processes.
- Validated for Gravity and Pre-Vac steam sterilization.

- Cleaning brush can be inserted through entire length of handle






e•Access™ Slide Lock Handle



+



Order Information

		Length	Catalog #
e•Access™ Slide Lock Handle (rotatable)		35 cm	ES8500SLR
		45 cm	ES8500-45SLR
Insert, Atraumatic Grasper (Wave)		35 cm	ES0537
		45 cm	ES0537-45
Insert, Atraumatic Grasper (Single Action)		35 cm	ES0570
		45 cm	ES0570-45
Insert, Traumatic Grasper		35 cm	ES0571
		45 cm	ES0571-45
Insert, Traumatic Grasper (Single Action)		35 cm	ES0573
		45 cm	ES0573-45
Insert, Fenestrated Bowel Grasper		35 cm	ES0574
		45 cm	ES0574-45